

thermostatic refrigerators

THERMOSTATIC REFRIGERATORS FOR BOD APPLICATIONS



Digital temperature control
via PID action with PT 100



External frame in steel sheet coated with epoxy powder.
Internal structure in ABS, with no joints.



Forced ventilation



Cooling system controlled
via an electronic valve



FTF

Technical data		FTF 180	FTF 220
Capacity	lt	180	260
Working temperature	°C	+4,0 ÷ +45,0	+4,0 ÷ +45,0
Resolution	°C	0,1	0,1
Precision	°C	±0,5	±0,5
Uniformity	°C	±1,0	±1,0
Heating power	W	130	450
Compressor	HP	1/5	1/4
Refrigerant Gas		R134a	R134a
Internal sockets 230V-50/60Hz	No.	2	4
Included Shelves	No.	3	4
Internal dimensions LxWxH	mm	430x500x680	480x500x1100
External dimensions LxWxH	mm	593x620x832	600x700x1358
Code		720.2072.18	720.2072.22

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 42