



LABORATORY
EQUIPMENT

melting points

WWW.FALCINSTRUMENTS.IT



MELTING POINTS

- 4 | Melting points for simple use
- 5 | Professional melting points

MELTING POINTS

Melting points for simple use



360 D



Designed for simultaneous measurement of three samples



Easy-to-use, internal LED light



Digital display: the two different rise speeds allow for quick measurement



Closed 100 capillaries D. 1-1.5 L 100 included

Technical data

360 D

Operating temperature	°C	Room ÷ +350
Precision at 20°C	°C	± 1,0
Precision at 300°C	°C	± 2,5
Display	n°	4 digit with 8 segments
Resolution	°C	0,1
Hold function		Yes
Number of samples	n°	3
Fast-speed ramp	°C/min	20 fixed
Slow-speed ramp	°C/min	from 1 to 10
Capacity	W	75
External dimensions LxWxH	mm	220x160x190
Weight	Kg	2,2
Code		621.0000.22

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 10

ACCESSORIES

Code	Capillaries Ø-mm	Length mm	Extremity	Pcs /Package
621.0000.28	1,0-1,5	80	opened	200
621.0000.30	1,0-1,5	100	opened	200
621.0000.38	1,0-1,5	80	closed	200
621.0000.40	1,0-1,5	100	closed	200

MELTING POINTS

Professional melting points



Designed for simultaneous measurement of three samples



Easy-to-use, internal LED light. The second display allows better visibility of the sample



Digital display: the two different rise speeds allow for quick measurement



Option to connect to an external printer. Forced cooling reduces analysis time. External manual control for easy and immediate HOLD function. Closed 100 capillaries D. 1-1.5 L 100 included



**360 D
COOLING**

Technical data

360 D COOLING

Operating temperature	°C	Room ÷ +350
Precision at 20°C	°C	± 1,0
Precision at 300°C	°C	± 2,5
Display	n°	4 digit with 8 segments
Resolution	°C	0,1
Hold function		Yes
Number of samples	n°	3
Fast-speed ramp	°C/min	25 fixed
Slow-speed ramp	°C/min	from 1 to 10
Capacity	W	75
External dimensions LxWxH	mm	220X160X190
Weight	Kg	2,2
Code		621.0000.23

V/HZ 220/240-50/60

SAFETY CLASS 0

PROTECTION CLASS IP 10

ACCESSORIES

Code	Capillaries Ø-mm	Length mm	Extremity	Pcs/Package
621.0000.28	1,0-1,5	80	open	200
621.0000.30	1,0-1,5	100	open	200
621.0000.38	1,0-1,5	80	closed	200
621.0000.40	1,0-1,5	100	closed	200



PRINTER

Code 621.0000.25

Notes

[illegible]

This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page, typical of notebook or legal stationery. There are no margins, text, or other markings present.



FALC INSTRUMENTS s.r.l.

Via G. M. Compagnoni, 2
24047 Treviglio (BG) - Italy
+39 0363 304660 | falc@falcinstruments.it

WWW.FALCINSTRUMENTS.IT