



MBO - 100 ml

These suspension mantles consist of a soft structure without metallic parts, while housings and interiors are in glasswire and glassfibre insulated.





Descrizione

These suspension mantles consist of a soft structure without metallic parts, while housings and interiors are in glasswire and glassfibre insulated.

Ideal for heating suspended heaters when traditional console versions cannot be used. A handy string allows heater fastening to mattress: simply roll it around the bowl.

MBO mantles are fitted with plug and 1.5 m. cord and with a 3 step power regulator with switch function.



Caratteristiche tecniche	
Power supply	V/HZ 230/50-60
Item	627.0624.01
Temperature	+300
Protection class	IP 0
Flask	ml 100
Bottom slot	no
Heating res. capacity	W 85