



MM - 1000 ml

MM heating mantles were developed by FALC on the base of long years of experience;





Descrizione

MM heating mantles were developed by FALC on the base of long years of experience;

- steel housings, oven coated with acidproof epoxy paint;
- semi - rigid structure that fits mattress shape drastically reducing risks of mattress breaking;
- thermal insulation guaranteed by ceramicfiber ;
- temperature control through a graduated thermoregulator.

All models are provided with ON /OFF switch with relevant green light.



Caratteristiche tecniche	
Power supply	V/HZ 230/50-60
Item	625.0500.10
Temperature	+350
External dimension (LxWxH)	mm 245x175
Weight	Kg 2.4
Protection class	IP 30
Flask	ml 1000
Precision on the plate	± 5
Heating res. capacity	W 450
Palloni Ø mm	131