



MA - 10000 ml

We attached great importance to safety by developing our MA mantles for large matrasses, provided with liquid recovery system, should matrass accidentally break.





Descrizione

We attached great importance to safety by developing our MA mantles for large mattresses, provided with liquid recovery system, should mattress accidentally break.

Housing and internal calottes completely made of stainless steel and covered with glassfibre cloth.

Resistances are accommodated in a ceramicfiber block, in order to make maximum use of the generated heat, while reducing outer wall temperature to a minimum. Temperature is controlled by a graduated thermoregulator.

MA mantles are fitted with an ON /OFF switch with relevant green light.



Caratteristiche tecniche	
Power supply	V/HZ 230/50-60
Item	625.0520.10
Temperature	+350
External dimension (LxWxH)	mm 460x330
Weight	Kg 11.0
Protection class	IP 44
Flask	ml 10.000
Precision on the plate	± 5
Heating res. capacity	W 1500
Palloni Ø mm	279