



Circulation cooling bath 15 lt - SBF 15T

The SBF 15 circulating bath with 15 liters tank is a device used to maintain a constant temperature of a liquid inside and outside a container. This type of bath is commonly used in laboratories to perform experiments requiring a constant temperature even below zero.

This instrument can be used in two ways:

- For thermostating closed external circuits from -30 to 150°C to power devices including photometers, refractometers, and viscometers or refrigeration coils, glass coolants.
- For thermostating applications in the bath from -30 to 150°C.

These circulating baths are employed in many applications in the chemical-pharmaceutical sector as well as in environmental analysis and material resistance testing.





Descrizione

- LCD touchscreen display
- Stainless steel 18/8 tank in a single mold
- Digital timer
- Eco-Friendly function, the cooling compressor operates only when necessary.
- Automatic compressor shutdown beyond 37.0°C
- Provision for simultaneous bath agitation and external circulation.
- Time starts counting when the instrument reaches the set temperature.
- Easy access for keeping the instrument clean and efficient.
- 4 language settings (Italian, English, French, Spanish)
- Program saving function (temperature/maintenance/pump speed)
- Programmable delayed start (day hour)



Caratteristiche tecniche	
Capacity	lt 15
Power supply	V/Hz 220-240/50-60
Item	612.1050.16
Temperature	from -30°C to + 150°C
Precision	± 0,2
Heating power	1300 W
Internal dimension (LxWxH)	mm 255x290x200
External dimension (LxWxH)	mm 300x440x730
Weight	Kg 41
Internal tank dimension	mm 250x190x200
Security class	DIN 12879 - 2
Protection class	IP 30
Refrigerant	W 550/600/845
Pump	premente
Pump pressure mbar 190	7,8 lt
Water head mt 1,8	Mt 2,5
Tank dimensions LxWxH	mm 250x190x200
Digital Timer	99:59:59 or ∞
Refrigerator power	-25°C_520 W / 0°C_600 W / 25°C_ 850W
Tub tap	No



Galleria immagini

