



AUTOCLAVE DA PAVIMENTO ATV 80 - 24 LT

The ATV 80 autoclave is a vertically loaded instrument for the sterilization of solids and liquids.

This model represents the perfect balance between compactness and efficiency thanks to its 24-liter chamber.

The basket is included.





Descrizione

- Sterilization cycles:
 - 4 pre-set cycles: 121°C - 135°C cycles for unwrapped solid and liquids
 - 1 customized cycle
- Pressure safety valve preset to 3 bar
- Over temperature protection
- Temperature monitoring
- Over current and short circuit protection
- Electric leakage protection
- Anti-scald chamber cover and bench
- Automatic troubleshooting
- Safety devices: mechanical lock to insert before starting the cycle
- Microcomputer control system monitors sterilization process.
- Saturated steam monitoring creates a sterilization environment and discharges the steam
- Built-in printer detecting the time, pressure and temperature parameters, serial number of the autoclave, the selected program, the progressive cycle counts, date, month, year
- Manual water filling



Caratteristiche tecniche	
Ambiente di lavoro	5-45°C, umidità relativa 10-85%
Pressione Di Lavoro	2.2 Mpa
Alimentazione	230 - Hz 50/60
Codice	608.1000.80
Dimensioni interne (LxPxA)	mm 243 x 500
Dimensioni esterne (LxPxA)	mm 500 x 450 x 950
Peso	Kg 60
Resistenza	W 2500
Volume camera	lt 24
Materiale camera	AISI 304 acciaio inox
Temperatura di sterilizzazione	100-136°C
Gamma di tempi di sterilizzazione	1-99 min
Potenza	2500 W
Acqua necessaria	demineralizzata o distillata
Dimensioni cestello	240x430 mm
Apertura coperchio	sollevare (su e giù)



Galleria immagini

