



Aluminum hotplat for industrial applications - XL-3

Ideal for heavy use (up to 40 kg) and material testing.
Large aluminium plate ensures uniformity of temperature.
Frame: the plate is not deformable.

Temperature is controlled by a thermostat with graduated scale in °C and expansion probe.





Descrizione

Ideal for heavy use (up to 40 kg) and material testing.
Large aluminium plate ensures uniformity of temperature.
Frame: the plate is not deformable.

Temperature is controlled by a thermostat with graduated scale in °C and expansion probe.



Caratteristiche tecniche

Temperature	from +50 to +350°C
Precision	+/-5°C
Heating power	2300W
External dimension (LxWxH)	500x470x200mm
Weight	17
Max stirring capacity	40 lt
Plate Diam.	400x500mm